

DrägerService®

Dräger

TECHNICAL SERVICE BULLETIN #168

Product: Narkomed 6000

Date: May 9, 2000

Re: Announcement of requirements to record software versions, CRC values and forward service reports to DrägerService.

Reference Doc:

Reason/Rational: To ensure proper software versions.

Solution/Procedure: All Narkomed 6000 software versions and CRC values shall be documented on the corresponding service report at the completion each repair. On a weekly or sooner basis all reports shall be forwarded electronically or a hard copy of the report shall be forwarded to DrägerService; Attention: Narkomed 6000 Technical Product Manager. Software and CRC configurations can be located on the About tab in the Setup notebook. Podbase (PBS) versions and CRC's must be verified for each pod installed and can be located in Primary Service mode in the Service Config notebook on the Update page.

Setup Notebook - About Page

	Version	CRC
NM6000	xxx	xxx
VPO	xxx	xxx
GAP	xxx	xxx
CV	xxx	xxx
SCR	xxx	xxx

Service Config Notebook - Update Page (PBS)

	Version	CRC
VPO	xxx	xxx
GAP	xxx	xxx
CV	xxx	xxx
SCR	xxx	xxx

Devices affected: All

Timing/When: All service calls.

Additional Info:

Cost:

Ordering Info/Parts:

Distribution: All Authorized Service Organizations, DrägerService

If you have any questions, please contact Technical Support by phone at 1-800-543-5047 or by fax at 1-215-721-5789.